

## Fitting Eyezen+ Lenses


Eyezen+ lenses require the following standard measurements:

- Monocular Distance PD
- Fitting Height (Simply measure from center pupil to the lowest point of the lens to acquire Fitting Height)

## VisionWeb™ Ordering Instructions

1. Enter **PD** and **Height**
2. Select **Single Vision** under the **Lens Type** drop-down menu
3. Select *Eyezen+* **Lens Design** (based on patient's age) and **Lens Material** desired
4. Select optional **Treatment**

**Bonus: Don't forget to add Crizal® No-Glare treatments to Eyezen+ Lenses for the most comprehensive protection.**



The screenshot shows the VisionWeb ordering interface with four numbered callouts:

- 1**: Points to the PD and Height input fields for both Right Eye (RE) and Left Eye (LE). The Dist. PD fields are set to 32.00 mm and the Height fields are set to 17.00 mm.
- 2**: Points to the **Lens Type** dropdown menu, which is set to **Single Vision**.
- 3**: Points to the **Lens Design** dropdown menu, which is set to **Essilor Eyezen+ 1 (Enter Height)**.
- 4**: Points to the **Available Treatments** list, which includes **Crizal Sapphire UV**.

Other visible elements include the **Job Type** dropdown set to **Frame To Come**, the **Lens Material** dropdown set to **Plastic 1.67**, and the **Available Treatments** list containing **Crizal Preventica**, **Crizal Sapphire UV**, **Crizal Avance UV**, **Crizal Alize UV**, and **Crizal Easy UV**.

Visit [EyezenPro.com](http://EyezenPro.com) to learn more.

Eyezen™+