






VARILUX®

LENSES AVAILABILITY for Eyecare Professionals

Effective **MAY 2018**

VARILUX® X DESIGN 		VARILUX® X FIT 		VARILUX® X 4D 		
MATERIAL	PRODUCT	INDEX	COLOR	POWER RANGE	ADD POWER*	MAX CYLINDER
PLASTIC	• 1.50 Index Clear	1.50		-10.00 to +6.00	0.75 to 4.00D	-6.00
	1.50 Index Essilor Polarized	1.50	Gray, Brown	-10.00 to +6.00	0.75 to 4.00D	-6.00
	1.50 Index Xperio UV™ Polarized	1.50	Gray, Brown	-10.00 to +6.00	0.75 to 4.00D	-6.00
	1.50 Index Transitions® Signature® VII	1.50	Gray, Brown	-10.00 to +6.00	0.75 to 4.00D	-6.00
	1.50 Index Transitions® XTRActive®	1.50	Gray, Brown	-10.00 to +6.00	0.75 to 4.00D	-6.00
TREXA® (TRIVEX®)	• Clear	1.53		-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® Signature® VII	1.53	Gray	-10.00 to +6.00	0.75 to 4.00D	-6.00
AIRWEAR® (POLYCARBONATE)	• Clear	1.59		-10.00 to +6.00	0.75 to 4.00D	-6.00
	Essilor Polarized	1.59	Gray, Brown	-10.00 to +6.00	0.75 to 4.00D	-6.00
	Xperio UV™ Polarized	1.59	Gray, Brown	-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® Signature® VII	1.59	Graphite Green, Gray, Brown	-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® XTRActive®	1.59	Graphite Green, Gray, Brown	-10.00 to +6.00	0.75 to 4.00D	-6.00
THIN&LITE®	• 1.67 Clear	1.67		-12.00 to +9.00	0.75 to 4.00D	-6.00
	1.67 Transitions® Signature® VII	1.67	Graphite Green, Gray, Brown	-12.00 to +9.00	0.75 to 4.00D	-6.00
	1.67 Transitions® XTRActive®	1.67	Graphite Green, Gray, Brown	-12.00 to +9.00	0.75 to 4.00D	-6.00
	• 1.74 Clear**	1.74		-14.00 to +13.00	0.75 to 4.00D	-6.00
	1.74 Transitions® Signature® VII	1.74	Gray, Brown	-14.00 to +13.00	0.75 to 4.00D	-6.00

Availability for Varilux S Series is the same as Varilux X Series.

MFH=14mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare coating is applied (for example, Crizal® No-Glare lenses).

*ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.50.

**Thin&Lite 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal Sapphire® 360°, Crizal Avancé UV™ or Crizal Alizé UV™.

Trivex is a registered trademark of PPG Industries Ohio, Inc.

- Available in Essential Blue Series™ lenses.



VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective **MAY 2018**

VARILUX® PHYSIO® W3+		TEMPORAL	NASAL	VARILUX® PHYSIO® W3+ FIT		TEMPORAL	NASAL	VARILUX® PHYSIO® W3+ EYECODE™		TEMPORAL	NASAL
MATERIAL	PRODUCT	INDEX	COLOR	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER				
PLASTIC	• 1.50 Index Clear	1.50		-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00				
	1.50 Essilor Polarized	1.50	Gray, Brown	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D 4.00D	-6.00				
	1.50 Xperio UV™ Polarized	1.50	Gray, Brown	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D 4.00D	-6.00				
	1.50 Index Transitions® Signature® VII	1.50	Gray, Brown	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00				
	1.50 Index Transitions® XTRActive®	1.50	Gray	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00				
TREXA® (TRIVEX®)	• Clear	1.53		-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00				
	Transitions® Signature® VII	1.53	Graphite Green, Gray, Brown	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00				
	Transitions® XTRActive®	1.53	Graphite Green, Gray, Brown	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00				
AIRWEAR® (POLYCARBONATE)	• Clear	1.59		-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00				
	Essilor Polarized	1.59	Gray, Brown	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D 4.00D	-6.00				
	Xperio UV™ Polarized	1.59	Gray, Brown	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D 4.00D	-6.00				
	Transitions® Signature® VII	1.59	Graphite Green, Gray, Brown	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00				
	Transitions® XTRActive®	1.59	Gray, Brown	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00				
THIN&LITE®	• 1.60 Clear	1.60		-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00D	-6.00				
	1.60 Transitions® Signature® VII	1.60	Gray, Brown	-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00D	-6.00				
	• 1.67 Clear	1.67		-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00				
	1.67 Essilor Polarized	1.67	Gray, Brown	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00				
	1.67 Xperio UV Polarized	1.67	Gray, Brown	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00				
	1.67 Transitions® Signature® VII	1.67	Gray, Brown	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00				
	1.67 Transitions® XTRActive®	1.67	Gray	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00				
	1.74 Clear**	1.74		-18.00 to +13.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00				
1.74 Transitions® Signature® VII	1.74	Gray	-18.00 to +13.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00					

MFH=14mm

VARILUX COMFORT® W2+		TEMPORAL	NASAL	VARILUX COMFORT® W2+ FIT		TEMPORAL	NASAL	VARILUX COMFORT® W2+ EYECODE™		TEMPORAL	NASAL
MATERIAL	PRODUCT	INDEX	COLOR	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER				
PLASTIC	• 1.50 Index Clear	1.50		-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00				
	1.50 Essilor Polarized	1.50	Gray, Brown	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00				
	1.50 Xperio UV Polarized	1.50	Gray, Brown	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00				
	1.50 Transitions® Signature® VII	1.50	Gray, Brown	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00				
	1.50 Transitions® XTRActive®	1.50	Gray	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00				
TREXA® (TRIVEX®)	• Clear	1.53		-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00				
	Transitions® Signature® VII	1.53	Gray, Brown	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00				
	Transitions® XTRActive®	1.53	Gray	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00				
AIRWEAR® (POLYCARBONATE)	• Clear	1.59		-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00				
	Essilor Polarized	1.59	Gray, Brown	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00				
	Xperio UV™ Polarized	1.59	Gray, Brown	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00				
	Transitions® Signature® VII	1.59	Graphite Green, Gray, Brown	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00				
	Transitions® XTRActive®	1.59	Gray, Brown	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00				
THIN&LITE®	• 1.60 Clear	1.60		-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00				
	1.60 Transitions® Signature® VII	1.60	Gray, Brown	-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00				
	• 1.67 Clear	1.67		-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00				
	1.67 Essilor Polarized	1.67	Gray, Brown	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00				
	1.67 Xperio UV Polarized	1.67	Gray, Brown	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00				
	1.67 Transitions® Signature® VII	1.67	Gray, Brown	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00				
	1.67 Transitions® XTRActive®	1.67	Gray	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.00D	-6.00				

MFH=14mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare coating is applied (for example, Crizal® No-Glare lenses).

*ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.00.

**Thin&Lite 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal Sapphire UV™, Crizal Avancé UV™ or Crizal Alizé UV™.

• Available in Essential Blue Series™ lenses.



VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective **MAY 2018**

VARILUX PHYSIO DRx™		LENS ID					
MATERIAL	PRODUCT	INDEX	COLOR	POWER RANGE	BASE CURVES	ADD POWER	CYLINDER
PLASTIC	• 1.50 Index Clear	1.50		-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50, 4.00D	-6.00
	1.50 Index Essilor Polarized	1.50	Gray, Brown, Gray-Green	-10.00 to +6.00	2.00, 3.25, 4.25, 6.25, 8.25	0.75 to 3.50D	-4.00
	1.50 Index Xperio UV™ Polarized	1.50	Gray, Brown, Gray-Green	-10.00 to +6.00	2.00, 3.25, 4.25, 6.25, 8.25	0.75 to 3.50D	-4.00
	1.50 Index Transitions® Signature® VII	1.50	Graphite Green, Gray, Brown	-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50, 4.00D	-6.00
	1.50 Index Transitions® XTRActive®	1.50	Graphite Green, Gray, Brown	-10.00 to +6.00	2.25, 4.25, 5.25, 8.25	0.75 to 3.50, 4.00D	-6.00
	1.50 Index Transitions® Vantage®	1.50	Gray	-10.00 to +6.00	2.25, 4.25, 5.25, 8.25	0.75 to 3.50, 4.00D	-6.00
TREXA® (TRIVEX®)	• Clear	1.53		-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® VII	1.53	Graphite Green, Gray, Brown	-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive®	1.53	Graphite Green, Gray, Brown	-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50, 4.00D	-6.00
AIRWEAR® (POLYCARBONATE)	• Clear	1.59		-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.59	Gray, Brown, Gray-Green	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.59	Gray, Brown, Gray-Green, Copper, Ruby	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® VII	1.59	Graphite Green, Gray, Brown	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive®	1.59	Graphite Green, Gray, Brown	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50, 4.00D	-6.00
	Transitions® Vantage®	1.59	Gray	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50, 4.00D	-6.00
THIN&LITE™	• 1.60 Clear	1.60		-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50, 4.00D	-6.00
	1.60 Transitions® Signature® VII	1.60	Gray, Brown	-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50, 4.00D	-6.00
	• 1.67 Clear	1.67		-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	1.67 Essilor Polarized	1.67	Gray, Brown, Gray-Green	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	1.67 Xperio UV™ Polarized	1.67	Gray, Brown, Gray-Green	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	1.67 Transitions® Signature® VII	1.67	Graphite Green, Gray, Brown	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	1.67 Transitions® XTRActive®	1.67	Graphite Green, Gray, Brown	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	• 1.74 Clear*	1.74		-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00
1.74 Transitions® Signature® VII	1.74	Gray	-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00	

MFH=17mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare coating is applied (for example, Crizal® No-Glare lenses).

*Thin&Lite 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal Sapphire UV™, Crizal Avancé UV™ or Crizal Alizé UV™.

Please see Xperio UV™ availability chart for full range of Xperio UV colors including Essilor® Sunwear Colors™.

- Available in Essential Blue Series™ lenses.

VARILUX PHYSIO SHORT DRx™



Rx availability for Varilux Physio Short DRx matches Varilux Physio DRx, except ADD power range is 0.75 to 3.00D across all materials.

MFH=14mm



VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective **MAY 2018**

VARILUX COMFORT DRx™		LENS ID						
MATERIAL	PRODUCT	INDEX	COLOR	POWER RANGE	BASE CURVES	ADD POWER	CYLINDER	
PLASTIC	• 1.50 Index Clear	1.50		-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50D	-6.00	
	1.50 Index Essilor Polarized	1.50	Gray, Brown, Gray-Green	-10.00 to +6.00	2.00, 3.25, 4.25, 6.25, 8.25	0.75 to 3.50D	-4.00	
	1.50 Index Xperio UV™ Polarized	1.50	Gray, Brown, Gray-Green	-10.00 to +6.00	2.00, 3.25, 4.25, 6.25, 8.25	0.75 to 3.50D	-4.00	
	1.50 Index Transitions® Signature® VII	1.50	Graphite Green, Gray, Brown	-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50D	-6.00	
	1.50 Index Transitions® XTRActive®	1.50	Graphite Green, Gray, Brown	-10.00 to +6.00	2.25, 4.25, 5.25, 8.25	0.75 to 3.50D	-6.00	
	1.50 Index Transitions® Vantage®	1.50	Gray	-10.00 to +6.00	2.25, 4.25, 5.25, 8.25	0.75 to 3.50D	-6.00	
TREXA® (TRIVEX®)	• Clear	1.53		-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50D	-6.00	
	Transitions® Signature® VII	1.53	Graphite Green, Gray, Brown	-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50D	-6.00	
	Transitions® XTRActive®	1.53	Graphite Green, Gray, Brown	-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50D	-6.00	
AIRWEAR® (POLYCARBONATE)	• Clear	1.59		-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50D	-6.00	
	Essilor Polarized	1.59	Gray, Brown, Gray-Green	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50D	-6.00	
	Xperio UV™ Polarized	1.59	Gray, Brown, Gray-Green, Copper, Ruby	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50D	-6.00	
	Transitions® Signature® VII	1.59	Graphite Green, Gray, Brown	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50D	-6.00	
	Transitions® XTRActive®	1.59	Graphite Green, Gray, Brown	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50D	-6.00	
	Transitions® Vantage®	1.59	Gray	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50D	-6.00	
THIN&LITE®	• 1.60 Clear	1.60		-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50D	-6.00	
	1.60 Transitions® Signature® VII	1.60	Gray, Brown	-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50D	-6.00	
	• 1.67 Clear	1.67		-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50D	-6.00	
	1.67 Essilor Polarized	1.67	Gray, Brown, Gray-Green	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50D	-6.00	
	1.67 Xperio UV™ Polarized	1.67	Gray, Brown, Gray-Green	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50D	-6.00	
	1.67 Transitions® Signature® VII	1.67	Graphite Green, Gray, Brown	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50D	-6.00	
	1.67 Transitions® XTRActive®	1.67	Graphite Green, Gray, Brown	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50D	-6.00	
	• 1.74 Clear*	1.74		-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50D	-6.00	
	1.74 Transitions® Signature® VII	1.74	Gray	-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50D	-6.00	

MFH=17mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare coating is applied (for example, Crizal® No-Glare lenses).

*Thin&Lite 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™ or Crizal Alizé UV™.

Please see Xperio UV™ availability chart for full range of Xperio UV colors including Essilor® Sunwear Colors™.

Trivex is a registered trademark of PPG Industries Ohio, Inc.

- Available in Essential Blue Series™ lenses.

VARILUX COMFORT SHORT DRx™

Rx availability for Varilux Comfort Short DRx matches Varilux Comfort DRx, except ADD power range is 0.75 to 3.00D across all materials.

MFH=14mm



Better Sight.
Better Life.

